

Website: www.eyetrusthq.ca



EyeTrust Equipment Price List 2022

Equipment	Features	Price*	ETN Price*
ET1 Lens Clock	 Features ball tips and measures +/- 0 to 20 diopters Packed in a protective box 	By Quote	20% off

Lens Clock

The lens clock is another critical measurement tool in the creation of prescription eyeglasses. The lens clock allows an optician, optometrist or ophthalmologist to measure the diopter curve of the spherical lens surface.

Features ball tips and measures +/- 0 to 20 diopters.

Packed in a protective box.





Item #: ET1



Website: www.eyetrusthq.ca

ET2 Cross Cylinder

Cross Cylinder sold individually

Measurements available: +/- 0.25, +/- 0.50, +/- 1.00

Price per cross cylinder

By Quote

20% off



eyeTrust 0

Cross Cylinder

3 Piece Set in a protective case Measurements: -0.25, -0.50, -1.00

Lenses are contrary dioptric in two axis positions and used to examine the degree and axis position of the cylindrical lens for determining astigmatism. To use, put the cross cylinder lens before cylindrical lens and make its one axis coincide with the axis position, then turn the cross cylinder lens 90° counter-clockwise and see the change of vision. If the patient's vision has no change, the degree of the cylindrical lens used can be considered as suitable, otherwise, the degree should be adjusted until patient can see correctly.



Item #: ET2



Website: www.eyetrusthq.ca

ET3 Lens Thickness Caliper

Test range: 0 - 20 mm Precision: 0.05 mm

Depth: 85 mm

eyeTrust 0

Packed in a protective black case

By Quote

20% off

Lens Thickness Caliper

The lens thickness caliper is a handheld, portable, and high precision measurement tool for determining lens thickness.

- 1. Test range: 0 20 mm
- 2. Precision: 0.05 mm
- 3. Depth: 85 mm

Packed in a protective black case.



Item #: ET3



Tel: 416 616 5192 Website: www.eyetrusthq.ca



ET4 Loose Prism Set 22 pieces

- This set includes 23 pieces that is stored in a durable aluminum case.
- 22 acrylic prism pieces (Size: 37 x 37mm):
- Diopter (Δ): 1/2, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, 20, 25, 30, 35, 40, 45 and 50
- 1 red filter lens

By Quote

20% off

Loose Prism Set 22 pieces

These high quality prisms are used to measure the deviation in strabismus. It is important to conduct polarized testing through prisms that are known to cause the least interference with polarization.

This set includes 23 pieces that is stored in a durable aluminum case.

22 acrylic prism pieces (Size: 37 x 37mm):

Diopter (Δ): 1/2, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, 20, 25, 30, 35, 40, 45 and 50

1 red filter lens



Item #: ET4



Website: www.eyetrusthq.ca

	- 4	//IG
eve	Frust	
0/0	1 450	

This set includes 17 pieces that is stored in a durable aluminum case. ET5 Loose Prism Set

16 acrylic prism pieces (Size: 37 x 37mm):

Diopter (Δ): 1/2, 1, 2, 4, 6, 8, 10, 12, 15, 20, 25, 30, 35, 40, 45 and 50

• 1 red filter lens

By Quote

20% off

Loose Prism Set 16 pieces

These high quality prisms are used to measure the deviation in strabismus. It is important to conduct polarized testing through prisms that are known to cause the least interference with polarization.

This set includes 17 pieces that is stored in a durable aluminum case.

16 acrylic prism pieces (Size: 37 x 37mm):

Diopter (Δ): 1/2, 1, 2, 4, 6, 8, 10, 12, 15, 20, 25, 30, 35, 40, 45 and 50

1 red filter lens

16 pieces



Item #: ET5



Website: www.eyetrusthq.ca

ET6 Prism bar set Horizontal and Vertical

- This set includes 2 prism bars that is stored in a durable aluminum case.
- 16 acrylic prism pieces (Size: 37 x 37mm):
- Diopter (Δ): 1/2, 1, 2, 4, 6, 8, 10, 12, 15, 20, 25, 30, 35, 40, 45 and 50

By Quote

20% off

Horizontal and Vertical Prism bar set


These high quality prisms are used to measure the deviation in strabismus.

This set includes 2 prism bars that is stored in a durable aluminum case.

16 acrylic prism pieces (Size: 37 x 37mm):

eyeTrust

Diopter (Δ): 1/2, 1, 2, 4, 6, 8, 10, 12, 15, 20, 25, 30, 35, 40, 45 and 50



Item #: ET6



Website: www.eyetrusthq.ca

ET7 Retinoscopy Rack

- Comes complete in its own foldable carrying case, which fits easily in pockets.
- Two very nice & durable retinoscopy bars.
- One is red and one is black.
- Each has 16 powers of lenses from: 0.5D-4D in half steps, 5D-10D in one diopter steps, 12D & 15D
- 10" x 1.75"

eyeTrust

By Quote

20% off

Retinoscopy Rack



Comes complete in its own foldable carrying case, which fits easily in pockets.

Two very nice & durable retinoscopy bars. One is red and one is black. Each has 16 powers of lenses from: 0.5D-4D in half steps 5D-10D in one diopter steps 12D & 15D

10" x 1.75"





Item #: ET7



Website: www.eyetrusthq.ca

ET8 Trial Lens Frame	 Range of PD scale: 50 - 80mm Single side (right or left) scale: 25 - 40mm Minimum scale Value: 1mm Axis scale step 5° Lens holders on right/left side available for 5 pieces of standard trial lenses inserted simultaneously Range of nose bridge adjustment: Height 0-23mm Angle: 0° - 360° Weight: 61g 	By Quote	20% off
----------------------	--	----------	---------

Trial Lens Frame



eyeTrust 0

Range of PD scale: 50 - 80mm Single side (right or left) scale: 25 - 40mm Minimum scale Value: 1mm

Axis scale step 5°

Lens holders on right/left side available for 5 pieces of standard trial lenses inserted simultaneously

Range of nose bridge adjustment: Height 0-23mm Angle: 0° - 360° Weight: 61g



Item #: ET8



Website: www.eyetrusthq.ca

 The Maddox Rod test is a specific test procedure to ascertain the amount of heterophoria or tropia, which is the ocular deviation of the allignment of the eyes. For large amounts of exophoria or esophoria, symptoms such as eyestrain, discomfort, headaches, browache, even double vision can occur. This product comes in a protective box with a white and red Maddox lens. Each lens has a silver handle that is 68mm in length. 	By Quote	20% off
---	----------	---------

Red & White Maddox Test

eyeTrust

The Maddox Rod test is a specific test procedure to ascertain the amount of heterophoria or tropia, which is the ocular deviation of the allignment of the eyes.

For large amounts of exophoria or esophoria, symptoms such as eyestrain, discomfort, headaches, browache, even double vision can occur.

This product comes in a protective box with a white and red Maddox lens. Each lens has a silver handle that is 68mm in length.



Item #: ET10



Website: www.eyetrusthq.ca

ET11C Optometry Tool Kit Set	 Comes complete in its own foldable carrying case 5 Pliers: Nose pad adjusting, straight needle nose, metal angling, dual nylon gripping, temple adjusting 4 Screwdrivers, 1 screw tweezer 1 Hex nut driver 2 Files: Flat and half round 1 Lens Caliper 1 Pen Light & 1 PD Ruler 1 Side Cutter 	By Quote	20% off
---------------------------------	--	----------	---------

Optometry Tool Kit Set V2

Comes complete in its own foldable carrying case and includes:

5 Pliers: Nose pad adjusting, straight needle nose, metal angling, dual nylon gripping, temple adjusting

4 Screwdrivers, 1 screw tweezer

eyeTrust 0

1 Hex nut driver

2 Files: Flat and half round

1 Lens Caliper

1 Pen Light & 1 PD Ruler

1 Side Cutter



Item #: ET11C



Website: www.eyetrusthq.ca

 Comes complete in its own foldable carrying case 4 Pliers: Nose pad adjusting, straight needle nose, metal angling, dual nylon gripping 4 Screwdrivers 1 Lens Caliper 1 Lens Clock 1 Pen Light 1 PD Ruler 	By Quote	20% off
---	----------	---------

Opticianry Tool Kit Set

eyeTrust

Comes complete in its own foldable carrying case and includes:

- 4 Pliers: Nose pad adjusting, straight needle nose, metal angling, dual nylon gripping
- 4 Screwdrivers
- 1 Lens Caliper
- 1 Lens Clock
- 1 Pen Light
- 1 PD Ruler



Item #: ET11B



Website: www.eyetrusthq.ca

	269 Piece Set, 38mm lenses.		
	Durable metal case.		
	• 0.25D steps to 6D.		
	• 0.50D steps to 10D.		
	• 1.00D to 16,18 & 20D.		
	 +/- Cyl 0.25D steps to 4D, 0.50D steps to 6D. 		
ET12 Trial Lens Set	 Prisms: 0.50D to 10D in 0.50D steps. 	By Quote	20% off
	Accessory Lenses: 1 red, 1 green, 1 black, 2 plano,		
	• 2 maddox, 1 slit.		
	• 1 polariscope, 2 cross, 1 frosted, 2 pinhole.		
	Set includes two Jackson Cross Cylinders.		
	• (+/-0.25D & +/-0.50D).		
	1 Lens Cleaning Cloth		



eyeTrust 0

Trial Lens Set

269 Piece Set, 38mm lenses. Durable metal case. 0.25D steps to 6D. 0.50D steps to 10D. 1.00D to 16,18 & 20D. +/- Cyl 0.25D steps to 4D, 0.50D steps to 6D. Prisms: 0.50D to 10D in 0.50D steps. Accessory Lenses: 1 red, 1 green, 1 black, 2 plano, 2 maddox, 1 slit. 1 polariscope, 2 cross, 1 frosted, 2 pinhole. Set includes two Jackson Cross Cylinders. (+/-0.25D & +/-0.50D). 1 Lens Cleaning Cloth



Item #: ET12



Website: www.eyetrusthq.ca

ET13 Flipper

- Power is indicated on the framework
- High quality plastic framework

Available powers: ±0.25D,±0.50D,±0.75D,±1.50D,±2.00D By Quote

20% off

Flipper (± 2.00D)


Power is indicated on the framework

High quality plastic framework

Available powers: ±0.25D,±0.50D,±0.75D,±1.50D,±2.00D

eyeTrust





Item #: ET13



Website: www.eyetrusthq.ca

Max Heating: 110° Minimum Heating: 50° Power of heating wires: 580W ET14 Frame Heater Power of electric fan: 6W 20% off By Quote ET-320 Noise: <50dB Wind speed: 1.8m/s Power supply: AC 110V/220V Weight: 1.5kgs

Frame Heater ET-320



eyeTrust

Max Heating: 110° Minimum Heating: 50°

Power of heating wires: 580W Power of electric fan: 6W

Noise: <50dB Wind speed: 1.8m/s

Power supply: AC 110V/220V

Weight: 1.5kgs







Item #: ET14

Website: www.eyetrusthq.ca

ET15 Frame Heater ET-340	 Metal box and brilliant painting endue a distinctive appearance A maximum heat airflow of up to 200 degrees celsius Temperature is invariable and has longer usage time Switch lights on to indicate on or off Built with a Nylon material black cover for extra durability and protection Voltage: AC220V/50Hz, AC110V/60Hz Temperature: Low - 170° High - 200° Single PTC heater: 130° - 150° Both PTC heater: 170° - 200° Power of PTC heater: 400W x2 Power of electric motor: 60W Temperature control: 2 selections of control Dimensions: L210 cm x D300 cm x H230 cm Weight: 3.2kgs 	By Quote	20% off

Frame Heater ET-340



eyeTrust 0

Voltage: AC220V/50Hz, AC110V/60Hz Temperature: Low - 170° High - 200°

Single PTC heater: 130° - 150° Both PTC heater: 170° - 200° Power of PTC heater: 400W x2 Power of electric motor: 60W

Temperature control: 2 selections of control Dimensions: L210 cm × D300 cm × H230 cm

Weight: 3.2kgs

. Metal box and brilliant painting endue a distinctive appearance

· A maximum heat airflow of up to 200 degrees celsius

· Temperature is invariable and has longer usage time

· Switch lights on to indicate on or off

. Built with a Nylon material black cover for extra durability and protection

Item #: ET15



Website: www.eyetrusthq.ca

•	Consecutive measurement, direct point sampling
	and high accuracy

- Uses mechanical hairspring that is leveled at the reflecting point of cornea
- Uses high-definition line-shaped sensor, advanced electronic system, and digital display, which makes test results visible, legible, and accurate
- LED light source and low power consumption
- Has +2.00D for eyesight compensation, accuracy self-check function, and VD measurement function

ET16 Pupilometer

eyeTrust

Range of Measurement:

Binocular Pupilary Distance: 45 - 82 mm Left or right Pupilar Distance: 22.5 - 41 mm

Indication Error: 0.5 mm Rounding Error: 0.5 mm Distance of target: 30 cm

Power Source: DC 3V, 2 AA batteries Time for automatic shut off: 1 minute

Size & Weight: 233 (L) x 165 (W) x 70 mm (H), 0.68

By Quote

90% off

Pupil Meter Range of Measurement: Binocular Pupilary Distance: 45 - 82 mm Left or right Pupilar Distance: 22.5 - 41 mm Indication Error: 0.5 mm Rounding Error: 0.5 mm Distance of target: 30 cm Power Source: DC 3V, 2 AA batteries Time for automatic shut off: 1 minute Size & Weight: 233 (L) x 165 (W) x 70 mm (H), 0.68 kg

· Consecutive measurement, direct point sampling and high accuracy

. Uses mechanical hairspring that is leveled at the reflecting point of cornea Uses high-definition line-shaped sensor, advanced electronic system, and digital display, which makes test results visible, legible, and accurate

LED light source and low power consumption

· Has +2.00D for eyesight compensation, accuracy self-check function, and VD measurement function

Item #: ET16



Tel: 416 616 5192 Website: www.eyetrusthq.ca

ET17 Auto- Lensometer	 Newest generation with Hartmann Sensor Technology Advanced simultaneous measurement of 108 points provides a faster measurement with great accuracy and reliability 5.6" colourful touch screen (TFT LCD) making it easier and convenient Can accurately and quickly measure the progressive lens Can display the pupil diameter and measure pupil distance Equipped with UV meter and built-in thermosensitive printer Fast measuring speed, clear interface, convenient operation, accuracy and good repeatability Please see image for full details 	By Quote	20% off
--------------------------	---	----------	---------



Auto-Lensometer

Specification:

Sphere: 0 to ±25D (0.01, 0.12, 0.25 D steps) Cylinder: 0 to ±10D (0.01, 0.12, 0.25 D steps) Axis: 0° ~ 180° (1° steps) ADD: 0 to 10D (0.01, 0.12, 0.25 D steps) Prism: 0 to 15Δ (0.01 / 0.12 / 0.25Δ steps) Step: 0.01 D Pupil Distance: 45 - 90 mm (0.25 mm steps)

Lens Diameter: 20 - 108 mm UV check: Outer UV (Optional) Abbe number range: 30 to 60

Measuring Mode:

Single, progressive, automatic recognition contact lens (soft and hard)

Speed: 0.1 seconds Display: 5.6" TFT LCD Printer: Thermal Printer

Power Supply: AC 100~120V/200~220V 50/60Hz Dimensions: 222(D)*161(W)*372(H)mm

Weight: 5 kgs (net)

· Newest generation with Hartmann Sensor Technology

eyeTrust

- · Advanced simultaneous measurement of 108 points provides a faster measurement with great accuracy and reliability
- 5.6" colourful touch screen (TFT LCD) making it easier and convenient
- · Can accurately and quickly measure the progressive lens
- Can display the pupil diameter and measure pupil distance
- · Equipped with UV meter and built-in thermosensitive printer
- · Fast measuring speed, clear interface, convenient operation, accuracy and good repeatability

Item #: ET17



Website: www.evetrusthg.ca



Specification:

- Range of Diopter measurement: 0 ~ ±25 D
- Minimum scale value: 0.125 D at $0 \sim \pm 5 D$, 0.25 D at ±5 D ~ 25 D
- Cylindrical lens axis of astigmatism: $0 \sim 180^{\circ}$, minimum scale value: 1°

ET18 Manual Lensometer

Prism Diopter:

- Pattern A:0~5 Δ , B:0~20 Δ , minimum scale value: 1 Δ
- Prism basal angle: 0~180°, minimum scale value: 1
- Adjustable range of eyepiece: ±5 D
- Measured size of lens: φ16 ~ φ80 mm
- Dimensions & Weight: 275(L)x130(W)x455(H) / 5.6kgs
- Illuminating lamp: 110V/220V (15W)

By Quote

20% off

Manual Lensometer

Specification:

Range of Diopter measurement: 0 ~ ±25 D Minimum scale value: 0.125 D at 0 \sim ±5 D, 0.25 D at ±5 D \sim 25 D Cylindrical lens axis of astigmatism: 0 ~ 180°, minimum scale value: 1°

Pattern A:0~5Δ, B:0~20Δ, minimum scale value: 1Δ Prism basal angle: 0~180°, minimum scale value: 1° Adjustable range of eyepiece: ±5 D Measured size of lens: φ16 ~ φ80 mm Dimensions & Weight: 275(L)x130(W)x455(H) / 5.6kgs Illuminating lamp: 110V/220V (15W)



Item #: ET18



Website: www.eyetrusthq.ca

ET20 High Grade Trial Lens Set	 266 Piece Set, Durable metal case. 0.25D steps to 6D. 0.50D steps to 10D. 1.00D to 16,18 & 20D. +/- Cyl 0.25D steps to 4D, 0.50D steps to 6D. Prisms: 0.50D to 10D in 0.50D steps. Accessory Lenses: 1 red, 1 green, 1 plano, 2 slit, 1 maddox, 1 crossed line, 1 frosted, 1 occluder, 1 +/- 0.50 cross cyl, 2 pinhole. Aperture of lenses: 36.5 mm Set includes 3 Jackson Cross Cylinders. (+/-0.25D, +/-0.50D & +/-1.00D). 1 Lens Cleaning Cloth 	By Quote	20% off
-----------------------------------	---	----------	---------

High Grade Trial Lens Set



eyeTrust 0

266 Piece Set, Durable metal case. 0.25D steps to 6D. 0.50D steps to 10D. 1.00D to 16,18 & 20D. +/- Cyl 0.25D steps to 4D, 0.50D steps to 6D. Prisms: 0.50D to 10D in 0.50D steps. Accessory Lenses: 1 red, 1 green, 1 plano, 2 slit, 1 maddox, 1 crossed line, 1 frosted, 1 occluder, 1 +/- 0.50 cross cyl, 2 pinhole. Aperture of lenses: 36.5 mm Set includes 3 Jackson Cross Cylinders. (+/-0.25D, +/-0.50D & +/-1.00D). 1 Lens Cleaning Cloth



Item #: ET20



Website: www.evetrusthg.ca



Specification:

- Range of Diopter measurement: 0 ~ ±20 D
- Minimum scale value: 0.125 D at $0 \sim \pm 5 D$, 0.25 D at ±5 D ~ 25 D
- Cylindrical lens axis of astigmatism: 0 ~ 180°, minimum scale value: 1°

Prism Diopter:

ET24 High Grade

Manual Lensometer

- Pattern A:0~5 Δ , B:0~20 Δ , minimum scale value: 1 Δ
- Prism basal angle: 0~180°, minimum scale value: 1
- Adjustable range of eyepiece: ±5 D
- Measured size of lens: ϕ 16 ~ ϕ 80 mm
- Dimensions & Weight: 275(L)x130(W)x455(H) / 6.4kgs
- Illuminating lamp: 110V/220V (15W) / Colours: Black or Jade Grey

By Quote

20% off

High Quality Manual Lensometer



Specification:

Range of Diopter measurement: 0 ~ ±20 D Minimum scale value: 0.125 D at 0 ~ ±5 D, 0.25 D at ±5 D ~ 25 D Cylindrical lens axis of astigmatism: 0 ~ 180°, minimum scale value: 1°

Prism Diopter:

Pattern A:0~5Δ, B:0~20Δ, minimum scale value: 1Δ Prism basal angle: 0~180°, minimum scale value: 1° Adjustable range of eyepiece: ±5 D Measured size of lens: φ16 ~ φ80 mm Dimensions & Weight: 275(L)x130(W)x455(H) / 6.4kgs Illuminating lamp: 110V/220V (15W) / Colours: Black or Jade Grey



Item #: ET24



Website: www.evetrusthg.ca

ET21 Phoroptor	 Auxiliary cylinder lenses (-2.00 D) Convergence adjustment Unique butterfly design Complete visual function testing Accurate measurements Comfortable and easy to use High grade coating of all optical lenses Please see image for full details Multiple options available 	By Quote	20% off
----------------	---	----------	---------



eyeTrust

Phoropter

Sphere Power: -19.00 D to +16.75 D (Step: 0.25 D) Cylinder Power: 0 D to 6.00 D (Step: 0.25 D)
Cylinder Axis: 0 to 180° (Step: 5°)

Cross Cylinder: ± 0.25 D Prism: 0 to 20△

PD Adjustment: 50 to 75 mm

Forehead rest Adjustment: 16 mm Vertex Distance: 13.75 mm

Dimensions: 335 (L) x 320 (W) x 98 (H) mm

Weight: Approx. 5 kg

- Auxiliary cylinder lenses (-2.00 D)
- Convergence adjustment
- Unique butterfly design
- · Complete visual function testing
- · Accurate measurements
- · Comfortable and easy to use
- High grade coating of all optical lenses

Item #: ET21



^{*}Additional cost for shipping and duties if applicable. All prices in CAD and tax (hst) applied.

^{*}eyeTrust Network (ETN) members get special discounted pricing

Equipment	Features	Price*	lmage
Occluder	Ask us for specific details	By Quote	On Request
Pinhole Occluder	Ask us for specific details	By Quote	On Request
Multi-Pinhole Occluder	Ask us for specific details	By Quote	On Request
PD Ruler & Occluder	Ask us for specific details	By Quote	On Request
PD Ruler	Black hard plastic	By Quote	On Request



Website: www.eyetrusthq.ca

	White measurement markings		
Screw Extractor	Ask us for specific details	By Quote	On Request
Stain Gauge	Ask us for specific details	By Quote	On Request
Portable Slit Lamp	Ask us for specific details	By Quote	On Request
20D Lens	Ask us for specific details	By Quote	On Request
90D Lens	Ask us for specific details	By Quote	On Request
Cross Cylinder (+/- 0.25)	Ask us for specific details	By Quote	On Request
Cross Cylinder (+/- 0.50)	Ask us for specific details	By Quote	On Request
Cross Cylinder (+/- 1.00)	Ask us for specific details	By Quote	On Request
Tonometer (Goldman)	Ask us for specific details	By Quote	On Request
Ultrasonic Cleaner	Ask us for specific detailsDifferent models available	By Quote	On Request
Vision Chart Screen	Digital screen – DS4236Ask us for specific details	By Quote	On Request
Vision Chart Screen	Digital screen – DS4901Ask us for specific details	By Quote	On Request
Hand Edge Grinder	Ask us for specific details	By Quote	On Request
Slit Lamp	Ask us for specific details	By Quote	On Request
BIO Bulbs	Ask us for specific details	By Quote	On Request
Topcon Slit Lamp bulb	Ask us for specific details	By Quote	On Request
Motorized Tables	Ask us for specific detailsMultiple models available	By Quote	On Request

^{*}Additional cost for shipping and duties if applicable. All prices in CAD and tax (hst) applied.

EyeTrust Optical Supply Price List 2022

eyeTrust 0

Optical Supply	Features	Price*	Image
Eyeglass Case	 Leather hard case (can select leather colour) Durable hinge Inside fabric embedded within the case (not glued like other case designs) Embossed logo on front outside of case Foil stamped logo on inside of case (can select foil stamp colour) 	Custom Quote	On Request Samples Available

^{*}eyeTrust Network (ETN) members get special discounted pricing



Website: www.eyetrusthq.ca

Cleaning Cloth	 Size: 16 x 16 or 16 x 19 cm Can choose colour Can choose edge pattern Can choose corner style Logo sample proof, design edits and printing included 200g, Microfiber cloth 	Custom Quote	On Request Samples Available
Lens Cleaner	 30, 60, 120 ml bottle available Alcohol free option available Solution colour optional Durable plastic bottle Spray head with cover 	Custom Quote	On Request Samples Available

^{*}Additional cost for shipping and duties if applicable. All prices in CAD and tax (hst) applied.

EyeTrust Stock Lens Price List 2022

eyeTrust 0

Lens Supply	Features	Price per pair*	Image
1.56 Stock Lenses	 1.56 CR39 Basic AR Basic scratch resistance coating 65 blank for + lenses (+0.25 to +6.00 / -0.25 to -2.00 CYL) 73 blank for - lenses (PL to -6.00 / -0.25 to -2.00 CYL) Ask us for specific details Add \$2 for -2.25-4.00 D CYL 	\$4	On Request Samples Available
1.59 Stock Lenses	 1.59 Polycarbonate Basic AR Basic scratch resistance coating 65 blank for + lenses (+0.25 to +6.00 / -0.25 to -2.00 CYL) 73 blank for - lenses (PL to -6.00 / -0.25 to -2.00 CYL) Add \$2 for -2.25-4.00 D CYL Ask us for specific details 	\$6	On Request Samples Available
1.59 Stock Lenses	 1.59 Polycarbonate No AR – Tintable Basic scratch resistance coating 	\$4	On Request

^{*}eyeTrust Network (ETN) members get special discounted pricing



Website: www.eyetrusthq.ca

	 65 blank for + lenses (+0.25 to +6.00 / -0.25 to -2.00 CYL) 73 blank for - lenses (PL to -6.00 / -0.25 to -2.00 CYL) Add \$2 for -2.25-4.00 D CYL Ask us for specific details 		Samples Available
1.59 Stock Lenses	 1.59 Polycarbonate Super AR Basic scratch resistance coating 65 blank for + lenses (+0.25 to +6.00 / -0.25 to -2.00 CYL) 73 blank for - lenses (PL to -6.00 / -0.25 to -2.00 CYL) Add \$2 for -2.25-4.00 D CYL Ask us for specific details 	\$8	On Request Samples Available
1.60 Stock Lenses	 1.60 Hi-Index Basic AR Basic scratch resistance coating 65 blank for + lenses (+0.25 to +6.00 / -0.25 to -2.00 CYL) 75 blank for - lenses (PL to -8.00 / -0.25 to -2.00 CYL) Add \$2 for -2.25-4.00 D CYL (only available from Pl to -6.00) Ask us for specific details 	\$10	On Request Samples Available

^{*}Additional cost for shipping and duties if applicable. All prices in CAD and tax (hst) applied.

eyeTrust 0

^{*}Min. order quantity for lenses is 2,000 in mixed order

^{*}eyeTrust Network (ETN) members get special discounted pricing